

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# GPA33 Fusion Protein



## Basic Information

**Catalog Number:**

Ag12655

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

TPQDVLRASQGKSVTLPCITYHTSTSSREGLIQWDKLLL  
THTERVVIWPF5NKNYIHGELYKNRVSISNNAEQSDA  
SITIDQLTMADNGTYECSVSLMSDLEGNTKSRVRLVL  
VPPSKPECGIEGETIIGNNIQLTCQSKEGSPTPQYSWK  
RYNILNQEQLAQPASGQPVSLKNI5TDTSGYYICTSS  
NEEGTQFCNITVAVRSPSMNVALYVGIAVGVVAALII  
IGIIYCC5CRGKDDNTE5KEDARPNREAYE5PPEQLR  
EL5RERE5DDYRQEEQR5TGRES5PDHLDQ  
(26-319aa encoded by BC074830)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.