

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PLAC4 protein



Basic Information

Catalog Number:

Ag12809

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MGTTTRSI PRD S P L G C L L S Y W S K F K F G L K K K K L I L Y C N T I
W V K Y K L D D E Y I W P K H G S I Y Y S D I L R L G L F C K R E G K W E E
I P Y V Q A F M A L Y W I T L L P G I R M P C I R D P H L A A P H R K F I S L
P R V S S V P Q S

(1-125 aa encoded by BC101615)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.