

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human IL-20 protein



Basic Information

Catalog Number:

Ag12828

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KTNLGSCVIATNLQEI R N G F S E I R G S V Q A K D G N I D I R I
L R R T E S L Q D T K P A N R C L L R H L L R L Y L D R V F K N Y Q T P D
H Y T L R K I S S L A N S F L T I K K D L R L C H A H M T C H C G E E A M K
K Y S Q I L S H F E K L E P Q A A V V K A L G E L D I L L Q W M E E T E
(26-176 aa encoded by BC074949)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.