

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

UBE4B Fusion Protein



Basic Information

Catalog Number:

Ag12836

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

EMIENHPLSTKLLVPSLMKFYTDVEHTGATSEFYDKFTI
RYHISTIFKSLWQNIHHGTFMEEFN SGKQFVRYINM
LINDTFFLLDESLESLKRIHEVQEEMKNKEQWDQLPRD
QQQARQSQAQDERVSRSYLALATETVDMFHILTQQ
VQKPFRLPELGPRLAAMLNFNLCGPKCRDLKVEN
PEKYGFEPKLLDQLTDIYQLDCARFAKAIADDQRSY
SKELFEEVISKMRKAGIKSTIAIEKFKLLAEKV EIVAKN
ARAEIDYSDAPDEFRDPLMDTLMTPVRLPSGTIMDR
SII LRHLLNSPTDPFNRQTLTESMLEPVPELKEQIQAW
MREKQNSDH

(825-1173 aa encoded by BC093696)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.