

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human NOP56 protein



Basic Information

Catalog Number:

Ag12878

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

DASRSSMGMDISAIIDLINIESFSSRVVLSSEYRQSLHTY
LRSKMSQVAPSLSALI GEAVGARLIAHAGSLTNLAKYP
ASTVQILGAEKALFRALKTRGNTPKYGLIFHSTFIGRAA
AKNKGRISRYLANKCSIASRIDCFSEVPTSVFGEKLR
QVEERLSFYETGEIPRKNLDVMKEAMVQAEAAAAEIT
RKLEKQEKRLKKEKRLAALALASSESSSTPEECEEM
SEKPKKKKQKPQEVQENGMEDPSISFSKPKKKKSF5
KEELMSSDLEETAGSTSIPKRRKSTPKEETVNDPEEAG
HRSGSKKKRKF5KEEPPVSSGPEEAVGKSSSKKKKFKHK
ASQED

(246-594aa encoded by BC104791)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.