

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RAD18 Fusion Protein



Basic Information

Catalog Number:

Ag12944

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MDSLAEsrwppGLAVMKTIDLLRCGICFEYFNIAmII
PQCshnyCSLCIRKFLSYKTQCPTCCVTVTEPDLKNN
RILDELVKSLNFARNHLLQFALESPAKSPASSSSKNLAV
KVYTPVASRQSLKQGSRLMDNFLIREMSGSTSELLIKE
NKSKFSPQKEASPAAKTKETRsvEEIAPDPSEAKRPEP
PSTSTLKQVTKVDCPVCgvNIPESHINKHLDSCLSREE
KKEsLRSSVHKRkPLPKTVYNLLSDRDlKkLkEHGLSI
QGnkQQLIKRHQEFVHMynaQCdALHPKsAAEIVQ
EIEIEKTRMRLEASKLNESVMVFTKDQTEKEIDEIHSK
YRkKHkSEFQLLVDQARKGYKkiAGMSQKTVTITKED
ESTEKLSsVCMGQEDNMTsvTNHFSQSKLDSPEELeP
DREEDSSsCIDIQEVLSSESdSCNssSSDIIRDLEEE
EAWEASHKNDLQDTEISPRQNRRTRAAESAEIEPRNK
RNRN

(1-495 aa encoded by BC001302)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.