

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human RAD18 protein



## Basic Information

**Catalog Number:**

Ag12954

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MDSLAESRWPPGLAVMKTIDDLLRCGICFEYFNIAIII  
PQC SHNYCSLCIRKFLSYKTQCPTCCVTVTEPDLKNN  
RILDELVKSLNFARNHLLQFALESPA KSPASSSSKNLAV  
KVYTPVASRQSLKQGSRLMDNFLIREMSGSTSELLIKE  
NKS KFSPQKEASPA AKTKETR SVEE IAPDPSEAKRPEP  
PSTSTLKQVTIKVDCPVC GVNIPESHINKHLD SCLSREE  
KKE SLRSSVHKRKPLPKTVYNLLSDRDLK K LKEHGLSI  
QG NKQQLIKRHQEFVHMYNAQC DALHPKSA AEIVQ  
EIENIEKTRMRLEASKL NESVMVFTKDQTEKEIDEIHSK  
YRKKHKSEFQLLV DQARKGYK KIAGMSQKT VTI TIKED  
ESTEKLSSVCMGQEDNMTSVTNHFSQSKLDSPEE LEP  
DREEDSSSIDIQEVLSSSES DSCNSSSSDIIRD LLEE  
EAWEASHKNDLQDTEISPRQNRRTRAAESAEIEPRNK  
RNRN

(1-495aa encoded by BC001302)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.