

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human NUP155 protein



## Basic Information

Catalog Number:

Ag12974

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, T-HIS, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

HGEPEEDIVGLQAFQERLNSYKCIDTLQELVNQSKA  
APQSPSPKPKPGPPVLSSDPNMLSNEEAGHHFEQML  
KLSQRSKDELFSIALYNWLIQVDLADKLLQVASPFLEP  
HLVRMAKVDQNRVRYMDLLWRYYEKNRSFSNAARVL  
SRLADMHSTEISLQQRLEYIARAILSAKSSTAISIAAD  
GEFLHELEEKMEVARIQLQIQETLQRQYSHSSVQDA  
VSQLDSELMDITKLYGEFADPFKLAECKLAIHCAGYS  
DPILVQTLWQDIIEKELSDSVTLSSSDRMHALSLKIVLL  
GKIYAGTPRFFPLDFIVQFLEQQVCTLNWDVGFVIQT  
MNEIGVPLRLLLEVYDQLFKSRDPFWNRMKKPLHLLDC  
IHVLLIRYVENPSQVLCERRRFTNLCLDAVCGYLVEL  
QSMSSSVAVQAITGNFKSLQAKLERLH  
(952-1391aa encoded by BC039257)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.