

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human GATA3 protein



Basic Information

Catalog Number:
Ag13089

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MEVTADQPRWVSHHHPAVLNGQHPDTHHPGLSHSY
MDAAQYPLPEEVDVLFNIDGQGNHVPPYYGNSVRAT
VQRYPPTHHGSQVCRPPLHGLPWLDGGKALGSHH
TASPWNLSPFKTSIHGSGPGLSVYPPASSSSLSGGH
ASPHLFTFPPTPPKDVSPDPSLSTPGSAGSARQDEKE
CLKYQVPLPDSMKLESSHSRGSMTALGGASSTHHPI
TTYPPYVPEYSSGLFPPSLLGGSPTGFGCKSRPKARSS
TEGRECVNCGATSTPLWRRDGTGHYLCNACGLYHKM
NGQNRPLIKPKRRLSAARRAGTSCANCQTTTTLWRR
NANGDPVCNACGLYYKLNINRPLTMKKEGIQTRNRK
MSSKSKKCKKVHDSLEDFPKNSSFNPAALSRHMSSLS
HISPFSSHMLTTPMHPPSSLSFGPHHPSSMVTA
M
(1-443 aa encoded by BC006793)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.