

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PDK1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag13121

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MRLARLLRGAALAGPGPLRAAGFSRSFSSDSSGSSPA  
SERGVPGQVDFYARFSPSPLSMKQFLDFGSVNACEK  
TSFMFLRQELPVRLANIMKEISLLPDNLLRTPSVQLVQS  
WYIQSLQELLDFKDKSAEDAKAIYDFDTVIRIRNRHN  
DVIPTMAQGVIEYKESFGVDPVTSQNVQYFLDRFYM  
SRISIRMLLNQHSLFLGGKGGSPSHRKHIGSINPNCN  
VLEVIKDGYENARRLCDLYIINSPELELEELNAKSPGQP  
IQVVYVPSHLYHMVFEFKNAMRATMEHHANRGVYP  
PIQVHVTLGNEDLTVKMSDRGGGVPLRKIDRLFNYMY  
STAPRPRVETSRAVPLAGFGYGLPISRLYAQYFQGDLLK  
LYSLEGYGTDAVIYIKALSTDSIERLPVYNKAAWKHYN  
(1-436 aa encoded by BC039158)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.