

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CD248 Fusion Protein



Basic Information

Catalog Number:

Ag13334

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

EPTQPPDFALAYRPSFPEDREPQIPYPEPTWPPPLSAP
RVVYHSSVLSVTRPVVVSATHPTLPSAHQPPVIPATH
PALS RDHQIPVIAANYPDLP S A Y Q P G I L S V S H S A Q P P
A H Q P P M I S T K Y P E L F P A H Q S P M F P D T R V A G T Q T T T H L
P G I P P N H A P L V T T L G A Q L P P Q A P D A L V L R T Q A T Q L P I I
P T A Q P S L T T S R S P V S P A H Q I S V P A A T Q P A A L P T L L P S
Q S P T N Q T S P I S P T H P H S K A P Q I P R E D G P S P K L A L W L P S
P A P T A A P T A L G E A G L A E H S Q R D D R W L
(399-687 aa encoded by BC051340)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.