

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human KLF8 protein



## Basic Information

Catalog Number:

Ag13365

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MVDMDKLI NNLEVQLNSEGGSMQVFKQVTSVRNR  
DPPEIEYRSNMTSPTLLDANPMENPALFNDIKIEPPEL  
LASDFSLPQVEPVDLDFHFKAPLQPASMLQAPIRPPK  
PQSSPQTLVVSTSTSDMSTSANIPTVLTPGSVLTSSQ  
STGSQQILHVIHTIPSVSLPNKMGGLKTI PVVVQSLPM  
VYTTLPADGGPAAITVPLIGGDGKNAGSVKVDPTSMS  
PLEIPSDSEESTIESGSSALQSLQLQEPAAMAQMQ  
GEESLDLKRRIHQCDFAGCSKVYTKSSHLKAHRIHT  
GEKPYKCTWDGCSWKFARSDELTRHFRKHTGIKPFRC  
TDCNRSFSRSDHLSLHRRRHDTM  
(1-359aa encoded by BC112109)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.