

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ZIC1 protein



## Basic Information

Catalog Number:

Ag13394

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

DFLFRNRGFGDAAAASAQHSLEAASAGGFGGPHG  
HTDAAAGHLLFPGLHEQAAGHASPNVVNGQMRGFS  
GDMYRPEQYGGVTSRSEHYAAPQLHGYPMNVN  
MAAHHGAGAFFRYMRQPIKQELICKWIEPEQLANPK  
KSCNKTFTMHVELVTHVTVHEVGGPEQSNHICFWEE  
CPREGKPFKAKYKLVNHIRVHTGEKPFPCPFPGCGKV  
FARSENKIHKRTHTGEKPFKCFEGCDRRFANSSDRK  
KHMHVHTSDKPYLCKMCDKSYTHPSSLRKHMKVHES  
SSQGSQSPAASSGYESSTPPTIVSPSTDNPTTSSLSP  
(99-447aa encoded by BC104848)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.