

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# CD36 Fusion Protein



## Basic Information

**Catalog Number:**

Ag13541

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MGCDRNCGLIAGAVIGAVLAVFGGILMPVGDLLIQKT  
IKKQVVLLEGTIAFKNWWKTGTEVYRQFWIFDVQNP  
QEVMMNSSNIQVKQRGPYTYRVFLAKENVTQDAE  
DNTVSFLQPNGAIFEPSLSVGTEADNFTVLNLAVAAA  
SHIYQNQFVQMLNSLHNKSKSMFQVRTLRELLWGY  
RDPFLSLVPYPVTTTVGLFYPYNNNTADGVYKVFNGKD  
NISKVAIIDTYKGRNLSYWESHCDMINGTDAASFPPF  
VEKSQVLQFFSSDICRSIYAVFESDVNLKIPVYRFVL  
PSKAFASPVENPDNYCFCTEKISKNCTSYGVLDISKC  
KEGRPVYISLPHFLYASPDVSEPIDGLNPNEEEHRTYLD  
IEPITGFTLQFAKRLQVNLVKPSEKIQLKLNLRNYIV  
PILWLNETGTIGDEKANMFRSQVTGKINLLGLIEMILLS  
VGVMFVAFMISYCACRSKTIK  
(1-472 aa encoded by BC008406)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.