

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human AGMAT protein



Basic Information

Catalog Number:

Ag13609

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PKG, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MLRLASGCARGPGVGVGARPAAGLFHPGRRQSRQA
SDAPRNQPPSPEFVARPVGVCSMMRLPVQTSPEGLD
AAFIGVPLDTGTSNRPGARFGPRRIEESVMLRTVNPS
TGALPFQSLMVADLGDVNVNLYNLQDSCRRIQEAYE
KIVAAGCIPLTGGDHTITYPIQAMAKKHGPGVGLLHV
DAHTDTTDKALGEKLYHGAPFRRCVDEGLDCKRVV
QIGIRGSSTLDPYRYNRSQGFRVVLAEDCWMKSLVP
LMGEVRQQMGGKPIYISFDIDALDPAYAPGTGTPEIA
GLTPSQ

(1-300 aa encoded by BC005090)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.