

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

PCGF6 Fusion Protein



Basic Information

Catalog Number:

Ag13611

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MEGVAVVTAGSVGAAKTEGAAALPPPPPPVSPPAL
TPAPAAAGEEGPAPLSETGAPGCSGSRPPELEPERSLGR
FRGRFEDEDEEEEEEEEEEEEEEDMSHFSRLEGG
RQDSEDEEERLINLSELPYILCSICKGYLIDATTITECL
HTFCKSCIVRHFYYSNRCPKCNIVVHQTQPLYNIRLDR
QLQDIVYKLVINLEEREKKQMHDFYKERGLEVPKPAVP
QPVPSSKGRSKKLVLESVFRIPPELDMSLLEFIGANEGT
GHFKPLEKFFVRVSGEATIGHVEKFLRRKMGLDPACQ
VDIICGDHLLSEQYQTLREIRRAIGDAAMQDGLLVLHYG
LVVSPLKIT

(1-352aa encoded by BC010235)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.