

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# BAG5 Fusion Protein



## Basic Information

Catalog Number:

Ag13701

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

FEFAQSLVREKIVPFYNGGNCVTDEFEEGIQDIIIRLT  
HVKTGGKISLRKARYHTLTKICAVQEIIEDCMKKQPSL  
PLSEDAHPSVAKINFVMCEV NKARGVLIALLMGVNNN  
ETCRHLSCVLSGLIADLDALDVCGRTEIRNYRREVVED  
INKLLKYLDLEEEADTTKAFDLRQNHSLKIEKVLKRMRE  
IKNELLQAQNPSELYLSSKTELQGLIGQLDEVSLKNPC  
IREARRRAVIEVQTLITYIDLKEALEKRKLFAC EEHPSHK  
AVWNVLGNLSEIQGEVLSFDGNRTDKNYIRLEELLTKQ  
LLALDAVDPOGEEKKAARKQAVRLAQNILSYLDLKS  
DEWEY

(139-488aa encoded by BC044216)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.