

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

APPL1 Fusion Protein



Basic Information

Catalog Number:

Ag13703

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MPGIDKLP~~IEETLEDSPQ~~TRSL~~LG~~VFEEDATAISNYMN
QLYQAMHRIYDAQNELSAATHLTSKLLKEYEKQRFPL
GGDDEVMSSTLQQFSKVIDELSSCHAVLSTQLADAM
MFPITQFKERDLKEILTLKEVFQIASNDHDAAINRYSRL
SKKRENDKVKYEVTEDVYTSRKKQHQTMMHYFCALN
TLQYK~~KKI~~ALLEPLLGYMQAQISFFKMGSENLEQLEE
FLANIGTSVQNVRRMDSDIETMQQTIEDLEVASDPL
YVPDPDPTKFPVNRNLRKAGYLNARNKTGLVSSTWD
RQFYFTQGGNLM~~SQ~~ARGDVAGGLAM
(1-323 aa encoded by BC028599)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.