

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human SNRNP200 protein



## Basic Information

Catalog Number:

Ag13741

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ELVEQTLSDLEQSKCISIEDEMDVAPLNLGMIAAYYI  
NYTTIELFSMSLNAKTKVRGLIEIISNAAEYENIPRHHE  
DNLLRQLAQKVPKLNPKFNDPHVKTNLLQAHLSR  
MQLSAELQSDTEEILSKAIRLIQACVDVLSNGWLSPA  
LAAMELAQMVTQAMWSKDSYLKQLPHFTSEHIKRCR  
DKGVESVFDIMEMEDEERNALLQLTDSQIADVARFC  
NRYPNIELSYEVVDKDSIRSGGPVVVLVQLEREEVETG  
PVIAPLFPQKREEGWWWVIGDAKSNLSIKRLTLQOK  
AKVKLDFVAPATGAHNYTLYFMSDAYMGCDQEYKFS  
VDVKEAETDSDSD

(C-350aa encoded by BC001417)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.