

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ABHD6 protein



## Basic Information

Catalog Number:

Ag13745

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

LWPSALIRIYYWYWRRTLGMQVRVHHEDYQFCYSF  
RGRPGHKPSILMLHGFSAHKDMWLSVVKFLPKNLHLV  
CVDMPGHEGTRSSLDLSDGQVKRIHQFVECLKLN  
KKPFHLVGTSMGGQVAGVYAAAYPSDVSSLCLVCPA  
GLQYSTDNQFVQRLKELQGSAAVEKIPLIPSTPEEMSE  
MLQLCSYVRFKVPQQILQGLVDVRIPHNNFYRKLFLFI  
VSEKSRYSLHQNMDKIKVPTQIIWKGQDQVLDVSGA  
DMLAKSIANCQVELLENCGHSVVMERPRKTAKLII DFL  
ASVHNTDNNKLD

(C-309aa encoded by BC001698)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.