

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

hIST1 Fusion Protein



Basic Information

Catalog Number:

Ag13874

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MLGSGFKAERLRVNLRLVINRLKLEKKKTELAQKARKE
IADYLAAGKDERARIRVEHIIREDYLVEAMEILELYCDLL
LARFGLIQSMKELDSGLAESVSTLIWAAPRLQSEVAEL
KIVADQLCAKYSKEYGKLCRTNQTIVNDRMLHKLSV
EAPPKILVERYLIEIAKNYNVPYEPDSVVMAEAPPGVE
TDLIDVGFDDVKKGGPGRGGSGGFTAPVGGPDGTV
PMPMPMPMPMSANTPFSYPLPKGPVDDINADKNISS
AQIVGPGPKPEASAKLPSRPADNYDNFVLPPELPSVPD
TLPTASAGASTSASEDIDFDLSRRFEELKKKT
(1-335 aa encoded by BC000116)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.