

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

C20orf4 Fusion Protein



Basic Information

Catalog Number:

Ag14039

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MAAVQMDPELAKRFFEGATVVIILNMPKGTEFGIDYN
SWEVGPKFRGVKMIPPGIHFLHYSSVDKANPKEVGPR
MGFFLSLHQGLTVLRWSTLREEVDLSPAPESEVEAM
RANLQELDQFLGPYPYATLKKWISLTNFISEATVEKLQ
PENRQICAFSDVLPVLSMKHTKDRVGQNLPRCGIECK
SYQEGLARLPEMKPRAGTEIRFSELPQMFPEGATPAE
ITKHSMDLSYALETVLNKQFPSSQDVLGELQFAFVC
FLLGNVYEA FEHWKRLNLLCRSEAAAMMKHHTLYINLI
SILYHQLGEIPADFFVDIVSQDNFLTSTLQVFFSSACS
IAVDATLRKKAKEKFAHLTKKFRWDFAAEPEDCAPV
VVLPPEGIEMG

(1-384 aa encoded by BC001751)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.