

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

EXD2 Fusion Protein



Basic Information

Catalog Number:

Ag14090

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

DDHSSWRKVLKQCQGVVDIPFRSKGMSRLGEEVNGE
ATESQQKPRNKKSMDGMVPGNHQGRDPRKHKRKPL
GVGYSARKSPDYDNCFLHAPDGGQLCTCDRRKAQWY
LDKIGELVSEEPFVVKLRFEPAGRPEPGDYLMVKE
NLCVVCGKRDSYIRKNVIPHEYRKHFPKEMKDHNSHD
VLLCTSCHAI SNYYDNHLKQQLAKEFQAPIGSEEGLR
LLEDPERRQVRSGARALLNAESLPTQRKEELLQALREF
YNTDVVTEMLQEAASLETRISNENYVPHGLKVQC
HSQGGLRSLMQLESRRWQHFLDSMQPKHLPQQWSV
DHNHQKLLRKFGEDLPQLS

(147-496 aa encoded by BC001962)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.