

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant rat Angptl4 protein



## Basic Information

Catalog Number:

Ag14117

Size:

50 µg

Form:

Available lyophilized

Species:

rat

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MRCAPTAGAALVLC AATAGLLSAQGRPAQPEPPRFA  
SWDEMNL LAHG LLQLGHGLREHVERTRGQLGALERR  
MAACGNACQGP KGTDPKDRVPEGQAPETLQSLQTQ  
LKAQNSKIQQ LFQKVAQQQR YLSKQNLRIQNLSQI  
DLLTPTHLDNGVDKTSR GKRLPKMAQLI GLTPNATRLH  
RPPRDCQELFQEG ERHSGLFQIQPLGSPFFLVNCEMT  
SDGGWTVIQRR LNSVDFNQSW EAYKDFGDPQGE  
FWLGLEKMHSITGDRGSQLAVQLQDWDGNAKLLQFP  
IHLGGEDTAYS LQLTEPTANELGATNVPNGLSLPFST  
WDQDHLRGLNCAKSLSGWWFGTCSHSNLNGQ  
YFHSIPRQRQRK KGI FWKTKWGRYYPLQATTLIQP  
MEATAAS

(1-405 aa encoded by BC078944)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.