

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human AKR1CL2 protein



Basic Information

Catalog Number:

Ag14381

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MGDIPAVGLSSWKASPGKVTEAVKEAIDAGYRHFDCA
YFYHNEREVGAGIRCKIKEGAVRREDLFIAKLWCTC
HKKSLVETACRKSLKALKLNLYLDLYLHWPMGFKPPHP
EWIMSCSELSFCLSHPRVQDLPLDESNMVIPSDTDFL
DTWEAMEDLVITGLVKNI GVSFNHEQLERLLNKPGL
RFKPLTNQIECHPYLTQKNLISFCQSRDVSVTAYRPLG
GSCEGVDLIDNPVIKRIAKEHGKSPAQILIRFQIQRNVI
VIPGSITPSHIKENIQVDFELTQHDMDNILSLNRNLRL
AMFPM

(1-307 aa encoded by BC002862)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.