

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant ATMPK8 protein



Basic Information

Catalog Number:

Ag14398

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

e coli-derived, PET28a, with N-terminal 6^oHis.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ERVPA PNGETA EESSDVERRAAAVASTLESEEADN

GGGYSARNLMKSASISGSKCIGVQSKTDKEDTIAEEE

DNETVAELTDKVASLHNS

(499-589aa encoded by NM_101675)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.