

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human IIP45 protein



Basic Information

Catalog Number:

Ag14609

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, T-HIS, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MVEAEELAQRLRLNLELLRQLWVGQDAVRRSVARAA
SESSLESSSYNSETPSTPETSSTSLSTSCPRGRSSVW
GPPDACRGDLRDVARSGVASLPPANCQHQESLGRPR
PHSAPSLGTSSLRDPEPSGRLGDPGPEAQTSRSILAQ
QSKLSKPRVTFSEESAVPERSWRLRPYLGWDWIAGSLD
TSSITSQPEAFFSKLQEFRETNKEECICSHPEPQLPGL
RESSGSGVEEDHECVYCYRVNRRLPVPVDPGTPCRL
CRTPRDQQGPGTLAQPAHVRVSIPLSILEPPHRYHIHR
RKSFDA SDTLALPRHCLLGWDIFPPKSEKSSAPRNLDL
WSSVSAEAQHQLSGTSSPFHPASPMQMLPPTPTW
SVPQVPRPHVPRQKP

(1-388 aa encoded by BC008068)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.