

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human TNK2 protein



Basic Information

Catalog Number:

Ag1472

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MQPEEGTGWLLELLSEVQLQQYFLRLRDDLNVTRLSH
FEYVKNEDLEKIGMGRPGQRRLWEAVKRRKALCKRKS
WMSKVFSGKRLEAEFPPHHSQSTFRKTSAPGGPAG
EGPLQSLTCLIGEKDLRLEKLDGGSFVVRRGEWDAP
SGKTVSVAVKCLKPDVLSQPEAMDDFIREVNAMHSL
DHRNLIIRLYGVVLTTPMKMVTTELAPLGSLLDRLRKHQG
HFLLGTLTRYAVQVAEGMGYLESKRFIHRDLAARNLLL
ATRDLVKIGDFGLMRALPQNDDHYVMQEHKVPFA
WCAPE

(N-term-300aa encoded by BC008884)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.