

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

NUP50 Fusion Protein



Basic Information

Catalog Number:

Ag14760

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MASEVLKNRAIKKAKRRNVGFESDTGGAFKGFKGLV
VPSSGGGRFSGFGSGAGGKPLEGLSNGNNITSAPPFAS
AKAAADPKVAFGSLAANGPTTLVDKVSNPKTNGDSQ
QPSSGLASSKACVGNAYHKQLAALNCSVRDWIVKH
VNTNPLCDLTPIFKDYEKYLANIEQQHGNSGRNSESES
NKVAAETQSPSLFGSTKLQQUESTFLFHGNKTEDTPDK
KMEVASEKKTDPSSLGATSASFNFGKKVDSSVLGSL
SVPLTGFSPGNSSLFGKDTTQSKPVSSPFPTKPLEG
QAEGDSGECGGDEEENDEPPKVVVTEVKEEDAFYS
KKCKLFYKDNFEKKEGIGTLHLKPTANQKTQLLVRAD
TNLGNILLNVLIPPNMPCRTTGKNNVLIVCVPNPPIDEK
(1-440aa encoded by BC028125)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.