

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# P2RY13 Fusion Protein



## Basic Information

Catalog Number:

Ag14779

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

ARVPYTHSQTNNKDCRLQNQLFIKETTFLAATNI

CMDPLIYIFLCKKFKTEKLPQMQRKTTASSQENHSSQ

TDNITLG

(253-333 aa encoded by BC041116)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.