

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human BHLHB9 protein



Basic Information

Catalog Number:

Ag15266

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KPQPVYEINEKNRPKDWSEVTIWPNAVTPAVLGFR
SQAPSEASPPSYIVLASEENACSLPVATACRPSRNT
RSCSQPIPECRFSDPCIQTIDEIRRQIRIREVNGIKPF
ACPCKMECYMDSEEFKLVSLKSTTDPLIHKIARIAM
GVHNVHPFAQEFINEVGVVTLIESLLSFSPEMRKKTIV
ITLNPPSGDERQRKIELHVVKHMCKETMSFPLNSPGQQ
SGLKILGQLTDFVHHYIVANYFSELFHLLSSGNCKTR
NLVLKLLLNSENPTAARDMINMKALAAKLIFNQKEA
KANLVSGVAIFINIKEHIRKGSIVVVDHLSYNTLMAIFR
EVKEIETM

(C-350aa encoded by BC041409)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.