

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human KIAA1430 protein



## Basic Information

Catalog Number:

Ag15295

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MDQFGDILEGEVDHSFFDSDFEEGKKCETNSVFDKQ  
NDDPKERIDKDTKNVNSNTGMQTENYLTEKGNERN  
VKFPPEHPVENDVTQTVSSFSLPASSRSKCLCDVTTG  
LKHVSI PNRI PKIVKEGEDDYTDGEESDDGKKYHVK  
SKSAKPSTNVKKSIRKKYCKVSSSSSSLSSSSSGSGTD  
CLDAGSDSHLSDSSPSSKSSKHKVSGITLLSPKHKYKS  
GIKSTETQPSSTAPKCGHYPEESEDVTDVSPSTPDI  
SPLQSFELGIANDQKVKKKQENVSQEIYEDVLDLKN  
NSKYLKAAAKGKEKHEPDVSSKSSVLDSSLDHRHKQ  
KVLHDTMDLNHLL

(1-350aa encoded by BC030535)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.