

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PPP1CC Fusion Protein



## Basic Information

**Catalog Number:**

Ag1550

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MADLDKLNIDSIIQRLLVVRGSKPGKNVQLQENEIRGL  
CLKSREIFLSQPILLELEAPLKICGDIHGQYYDLLRLFY  
GGFPPESNYLFLGDYVDRGKQSLETICLLAYKIKYPEN  
FFLLRGNHECASINRIYGFYDECKRRYNIKWKFTDC  
FNCLPIAAIVDEKIFCCHGGLSPDLQSMEQIRRI MRPT  
DVDPQGLLCDLLWSDPKDVLVWGENDRGVSFTFG  
AEVVAKFLHKHDLDLICRAHQVVEDGYEFFAKRQLVT  
LFSAPNYCGEFDNAGAMMSVDETLMCSFQILKPAEK  
KKPNATRPVTPPRGMITKQAKK  
(1-323 aa encoded by BC014073)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.