

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

KIAA1199 Fusion Protein



Basic Information

Catalog Number:

Ag15527

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

TVAAGCPDQSPQLQPWNPGHDQDHHVHIGQGKTL
LTSSATVYSIHISEGGKLVKIDHDEPIVLRTRHILIDNGG
ELHAGSALCPFQGNFTIILYGRADEGIQDPYGLKYI
GVGKGALELHGQKLSWTFLNKTLHPGGMAEGGYF
FERSWGHGRGIVHVVDPKSGTVIHSRFDYRSKKESE
RLVQYLNVPDGRILSVAVNDEGSRNLDDMARKAMT
KLGSKHFLHLGFRHPWSFLT VKGNPSSSVEDHIEYHG
HRGSAARVFKLFQTEHGEYFNVLSSEWVQDVEW
TEWFDHDKVSQTKGGEKISDLWKAHPGKICNRPIDIQ
ATTMDGVNLSDEVVYKKGQDYRFA
(31-386 aa encoded by BC020256)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.