

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human TIGD7 protein



Basic Information

Catalog Number:

Ag15597

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

DI~~C~~LP~~G~~KK~~I~~NKERLSAFLCANADGTHKLSIIIGKSKLPK
SVKEDTSTLPVIYKPSKDVWFTRELFSEWFFQNFVPE
VRHFQLNVLRFHDEDVRA~~L~~LLDSCPAHPSESLS~~E~~D
GRIKCMFFPHNTSTLIQPMNQGVILSCKRLYRWKQLE
ESLVIFEESDDEQEKGDKGVSKIYNIKSAIFNWAKS
WEEVKQITIANAWENLLYKKEPEYDFQGLEHGDYREIL
EKCGELET~~K~~DDDRVWLN~~G~~DEEK~~G~~CLLTKGGITKEV
VQKGGAEKQTA~~E~~FKLSAVRESLDYLLDFV~~D~~ATPEFQ
RFHFTLKEMQQEIVKQ~~F~~QSKIHSRIGSFLKPRPHNIK
DSFSGPSTSGSNH

(C-353aa encoded by BC069623)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.