

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human NOSIP protein



Basic Information

Catalog Number:

Ag1562

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MTRHGKNCTAGAVVYTYHEKKKDTAASGYGTQNI
RSLRDAVKDFDCCCLSLQPC HDPVVTPDGYL
YERAILEYI LHQKKEIARQMKAYEKQ
RGRTRREEQKELQRAASQDH V
RGFLEKESAI VSRPLNPFTAKALS
GTSPPDDVQPGPSV GPPSKDKDK
VLPSPFWIPSLTPEAKATKLEKPSRT
VTCP MSGKPLRMSDLTPVHFTPLD
SSVDRVGLITRSERYVCA VTRD
SLSNATPCAVLRPSGAVVTLECV
EKLIRKDMVD PVTGDKLTDRDII
VLQRGGTGFA GSGVKLQAEKSRPV
MQA

(1-301 aa encoded by BC010077)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.