

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

OPLAH Fusion Protein



Basic Information

Catalog Number:

Ag15676

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

ELIGQYGLDVVQAYMGHIQANAELAVRDM LRAFGTS
RQARGLPLEVSS EDHMDDGSP IRLRVQISLSQGS AVF
DFSGTGPEVFGNLNAPRAVTL S ALIYCLRCLVGRDIPL
NQGC LAPVRVVI PRGSILDP SPEAAV VGGNVLT SQRV
VDVILGAFGACAASQ GCMNNVTLGNAHMGYYETVA
GGAGAGPSWHGRSGVHSHMTNTRITDPEILES RYPVI
LRRFELRRGSGGRGRFRGGDGV TRELLFREEALLSVL TE
RRAFRPYGLHGGE PGARGLNLLIRKNGRTVNLGGKTSV
TVYPGDVFC LHTPGGGGYGDPE DPAPPPGSP PQALA
FPEHGSVYEYRRAQEAV

(C-350aa encoded by BC142672)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.