

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

EBF3 Fusion Protein



Basic Information

Catalog Number:

Ag15693

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

CNQNCLKNAGNPRDMRRFQVVVSTTVNVDGHVLAV
SDNMFVHNNSKHGRRARRLDPSEATPCIKAI SPSEGW
TTGGATVIIIIGDNFFDGLQVVF GTMLVWSELITPHAIR
VQTPPRHIPGVVEVTL SYKSKQFC KGAPGRFVY TALN
EPTIDYGFQRLQKVI PRHPGDPERLPKEVLLKRAADLV
EALYGMPHNNQEII LKRAADIAEALYSVPRNHNQIPTL
GNNPAHTGMMGVNSFSSQLAVNVSETSQANDQVG
YSRNTSSVSPRGVVPSTPQQSNYNTVSTSMNGYGS
GAMASLGVPGSPGLNGSSANSPYGMKQKSAFAPVV
RPQASPPPSCSANGNGLQAM

(194-543 aa encoded by BC126130)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.