

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human IL-18RAP protein



Basic Information

Catalog Number:

Ag15927

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KPDILDPVEDTLEVELGKPLTISCKARFGFERVFNPIK
WYIKSDLEWEVSVPEAKSIKSTLKDEIIERNIILEKVT
QRDLRRKFVCFVQNSIGNTTQSVQLKEKRGVLLYILL
GTIGTLVAVLAASALLYRHWIEIVLLYRTYQSKDQTLG
DKKDFDAFVSYAKWSSFPSEATSSLSEEHLALSFPDV
LENKYGYSCLLERDVAPEGVYAEDIVSIIKRSRRGIFIL
SPNYVNGPSIFELQAAVNLAADDQTLKLIKFCYFQE
PESLPHLVKKALRVLPTVTWRGLKSVPPNSRFWAKMR
YHMPVKNSQGFTWNQLRITSRIFQWKGLSRTETTGRS
SQPKEW

(250-599 aa encoded by BC069630)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.