

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SYCP1 Fusion Protein



Basic Information

Catalog Number:

Ag16005

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

GTAESKQLNVYEIKVNKLELELES AKQKFG EITDTYQK
EIEDKKISEENLLEEVEKAKVIADEAVK LQKEIDKRCQH
KIAEMVALMEKHKHQYDKII EERDSELGLYSKEQE QS
SLRASLEIELSNLKA ELLSVKKQLEI EREEKEK LKREAKE
NTATLKEKDKKTQTF LLETPEIYWK LDSKAVPSQTVS
RNFTSVDHGISKDKRDY LWTSAKNTLSTPLPKAYTVKT
PTKPKLQ QRENLNPI EESK KKRKMAFEFDINS DSSETT
DLLSMVSEETLKTLYRNN NPPASHLCVKT PKKAPSSL
TTPGSTLKF GAIRKMREDR WAVIAKMDRKKK LKEAEKL
FV

(629-976 aa encoded by BC126266)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.