

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# C1orf107 Fusion Protein



## Basic Information

**Catalog Number:**

Ag16199

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

LKRPEYEA V F V G N I D D H F R I G V A I L Q R S I R L Y A P F Y S S  
D I L I A S P L G R T I I G G E G E K K R D F D L S S I E L L I D Q A D I Y L  
M Q N W E H V L H L M N H M N L L P L D S H G V D F S R V R M W S L N  
N W S K Y Y R Q T L L F G A L Q D A Q I N S V F N K Y C V N M Q G Q V  
A V R N V P M T G S I S H V L V Q L P H V F Q R M E A E N L A S V I D A  
R F N F F V N K I L P Q Y R D A V M S H T L I Y I P S Y F D F V R L R N Y F K  
K E E L N F T H I C E Y T Q K S G V S R A R H F F L Q G E K Q F L L F T E R  
F H F Y K R Y T I K G I R N L I F Y E L P T Y P H F Y S E I C N M L R A T N R G  
E E A T W T C T V L Y S K Y D A Q R L A A V V G V E R A A Q M L Q S N  
K N V H L F I T G E K

(C-350aa encoded by BC022964)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.