

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human SEMA5A protein



Basic Information

Catalog Number:

Ag16346

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

VCVFNLSAIAQAFSGPFKYQENSRSAWLPYPNP
FQC GTVDQGLYVNLTERNLQDAQKFLMHEVVQPV
TVPSFMEDNSRFSHVAVDVVQGREALVHIIYLA
TDYG TIKKVRVPLNQTSSSCLLEEIELFPERRR
PIRSLQILHS QSVLFVGLREHVVKIPLKRCQ
FYRTRSTCIGA QDPYC GWDVVMKKCTS
LEESLSMTQWEQSISACPTRNLTVD GHF
GVWSPWTPCTHTDGS AVGSLC RTRSCDSPAP
QCGGWQCEGPGMEIANCSRNGGWTPWTSW
SPCST TCIGIGFQVRQRSCSNPTPRHGGRV
CVGQNR EERYCN EHLLCPPHMFWTGWGP
WERCTAQ (319-668aa encoded by BC115696)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.