

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# CYP2A6 Fusion Protein



## Basic Information

Catalog Number:

Ag16377

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

FLSLLRMMLGIFQFTSTSTGQLYEMFSSVMKHLPGPQ  
QQAFQLLQGLDFIAKKVEHNQRTLDPNSPRDFIDSF  
LIRMQEEEEKNPTEFYLNLMVMTTLNLFIGGTETVSTTL  
RYGFLLLMKHPEVEAKVHEEIDRVIGKNRQPKFEDRAK  
MPYMEAVIHEIQRFQDVI PMSLARRVKKDTKFRDFFLP  
KGTVEVYPMGLSVLRDPSFFSNPQDFNPQHFLNEKGQ  
FKKSDAFVPPSIGKRNCFGEGLARMEFLFFTTVMQN  
FRLKSSQSPKIDVSPKHVGFATIPRNYTMSFLPR  
(198-494 aa encoded by BC096253)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.