

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human KIAA1383 protein



## Basic Information

Catalog Number:

Ag16447

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

PQLPEDKYLSDASFTENS DTSRQISGVFDEPSTSKET  
KLIKYTEKKTVDVDCSKNRINNVSL EEVVSPANSIIPERLT  
PTNILGGNVEMKIQSPCVFQQDAVVDRIVDK EIDIRQ  
VKTTDNDILMADISDKRTGKN SCYENISELKYSDDLSS  
PCYSEDFCTSEDTSRSFKAHDSSRTENPKHSQYTSK  
SSDTGVSKKKNSSDRSSILSPPFSA GSPVHSYRKFHISK  
TQDKSLEEASSISADLSSTHWTEQKENQIDQNSMH  
NSEITKRAQDISVKTRSSWKSLEK SQSPQTSQVSSYLP  
SNVSELNVLDSSTSDHFEEGNDDV GSLNISKQCKDIC  
ELVINKLPGYTM

(555-905 aa encoded by BC109122)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.