

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ARHGAP21 Fusion Protein



## Basic Information

**Catalog Number:**

Ag16515

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

LDSSRLSPEVQSVAESKGD EADDERSELI SEGRPVETD  
SESEFPVFPTALTSE RLF RGLQEVT KSSRRN SEGS ELS  
CTEGSLTSSLDSRRQLFSSHKLIECDTL SRKKSARFKSD  
SGSLGDAKNEKEAPSLTKVFDVMKKGKSTGSLLPTRG  
ESEKQ EPTWTKIADRLKLRPRAPDDMFGVGNHKV  
NAETA KRKSIRRRHTLGGHRDATEISVLNFWK VHEQS  
GERESEL SAVNRLKPKCSAQDLSISDWLARELRRTSTS  
DLSRGEIGDPQTENPSTREIATTDTPLSLHCNTGSSSS  
TLASTNRPLLSIPPQSPDQINGESFQNVSKNASSAAN  
AQPHKLETPGTKA EFHPCL

(1598-1957 aa encoded by BC118915)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.