

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ADAM8 protein



Basic Information

Catalog Number:

Ag16524

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

EGGEGGRHAVYQAEHLLQTAGTCGVSDDSLGSLLGP
RTAAVFRPRPGDSLPSRETRYVELYVVVDNAEFQMLG
SEAAVHRVLEVVNHVDKLYQKLNFRVVLVGLIWNLS
QDRFHVSPDPSVTLELLTWQARQRTRRLHDNVQL
ITGVDFGTTVGFARVSAMCSHSSGAVNQDHSKNP
VGVACTMAHEMGNLGMHDENVQGCRCQERFEA
GRCIMAGSIGSSFPRMFSDCSQAYLESFLERPQSVCL
ANAPDLSHLVGGPVCGNLFVERGEQCDCGPPEDCRN
RCCNSTTCQLAEGAQCAHGTCCQECKVKPAGELCR
PKKDMCDLEEFCDGRHPECPEDAFQEN
(145-494aa encoded by BC115404)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.