

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TIGD7 Fusion Protein



Basic Information

Catalog Number:

Ag16539

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

DI~~C~~LP~~G~~KK~~I~~NKERLSAFLCANADGTH~~H~~KLKSI~~I~~IGKSKLPK
SVKEDTSTLPVIYKPSKDVWFTRELFSEWFFQNFVPE
VRHFQLNVLRFHDEDVRA~~L~~LLLDSCPAHPSES~~L~~TSED
GRIKCMFFPHNTSTLIQPMNQGVILSCKRLYRWKQLE
ESLVIFEESDDEQEKGD~~K~~GVSKI~~K~~IYNIKSAIFNWAKS
WEEVKQITIANAWENLLYKKEPEYDFQGLEHGDYREIL
EKCGELET~~K~~DDDRVWLN~~G~~DEEK~~G~~CLLKTGGITKEV
VQKGGAEKQTAEFKLSAVRESLDYLLDFVDATPEFQ
RFHFTLKEMQQEIVKQ~~F~~QSKIHSRIGSFLKPRPHNIK
DSFSGPSTSGSNH

(197-549 aa encoded by BC069623)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.