

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human IL-17RD protein



## Basic Information

Catalog Number:

Ag16636

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ALDLWEDFSLCREGQREWWIQKHESQFIIVVCSKGM  
KYFVDKKNYKHKGGGRGSGKGLFLVAVSAIAEKLKQ  
AKQSSAALSFKIAVYFDYSCEGDVPGILDSTKYRLM  
DNLPLQLCSHLHSDHGLQEPGQHTROGSRNRYFRSK  
SGRSLYVAICNMHQFIDEPPDWFQFVFPFHPPPLRY  
REPVLEKFDSGLVLDNDVMCKPGPESDFCLKVEAAVLG  
ATGPADSQHESQHGGLDQDGEARPALDGSAAALQPLL  
HTVKA GSPSDMPRDSGIYDSSVPSSELSPLMEGLSTD  
QTETSSLTESVSSSSGLGEEPPALPSKLLSSGCKADL  
GCRSYDELHAVAPL

(358-707 aa encoded by BC111702)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.