

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human MED26 protein



Basic Information

Catalog Number:

Ag16915

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

QQLDRVDETPGPPHPKPPRCFSFSPRNSRHEGSFAR
QQSLYAPKGSVPPSPRPQALDATQVPSPLPLAQPST
PPVRRLELLPSAESPVCWLEQPESHQRLAGPGCKAGL
SPAEPLLSLAGFSPDSSKADSDAASSGGSDSKKKRYR
PRDYTVNLDGQVAEAGVKPVRLKERKLTFFDPMTRQIK
PLTQKEPVRA DSPVHMEQQSRTELDKQEAKASLQSP
FEQTNWKELSRNEIQSYLSRQSSLLSSGAQTPGAH
HFMSYLLKQEESTRQGARQLHVLVLPQSPPTDLPGLTR
EVTQDDLDRIQASQWPGVNGCQDTQGNWYDWTQ
CISLDPHGDDGRLNLPVYVCLD
(251-600 aa encoded by BC110430)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.