

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ZNF24 protein



Basic Information

Catalog Number:

Ag1724

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PKG, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSAQSVEEDSILIIPTPDEEEKILRVKLEEDPDGEEGSSI
PWNHLPDPEIFRQRFQFGYQDSPGPRAVSQLREL
RLWLRPETHTEQIILELVVLEQFVAILPKELQTWVRDH
HPENGEEAVTVLEDLESELDDPGQPVSLRRRKREVLVE
DMVSQEEAQGLPSELD AVENQLK WASWELHSLRH
CDDDGRTENGALAPKQELPSALESHEVPGLSMGVP
QIFKYGETCFPKGRFERKRNP SRKKQHICDECGKHFS
QGSALILHQRIHSGEKPYGCVCEGKA FSRSSILVQH
RVHTGEKPYKLECGKAFSQNSGLINHQRHTGEKPY
ECVQCCKSYSSNLF RHQRRHNAEKLNNVVKV
(1-368 aa encoded by BC016801)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.